

## Model No. 2535-II

This machine has been developed for rapidly disintegrating wet pulp (concentration 20 to 30%). Machine configuration: the material is charged into special agitating impellers that rotate at a high speed; then the pulp is disintegrated at an instant and discharged from the lower port. Both dry pulp and wet pulp can be processed, remaining less grains. Since the specimen tank can be removed, efficient measurement can be made.

Motor 100/110VAC 0.2kW steplessly variable

100/110VAC 50/60Hz 5A **Power Source Outer Dimensions** 270×530×1,160mm

**Instrument Weight** 60kg

