

**Pulp Disintegrator  
(with Stepless Speed Changer)**

**Model No. 2529**

This machine disintegrates the pulp into single fibers, by rotating the pulp with water at a high speed. This disintegration method is specified by JIS and ISO. Its applicable range is wide: paper sheet making, preparation of specimen for measurement of freeness, bleaching and screening. The rotation speed can be changed as desired. It is used for operation in a condition different from disintegration specified by the standard.

**Specifications**

<b>Capacity:</b>	2 liters
<b>Concentration:</b>	max. 3%
<b>Rotation Speed:</b>	750 to 3,000rpm
<b>Disintegration Time:</b>	set by the timer
<b>Optional:</b>	jacket type specimen tank
<b>Motor:</b>	three-phase 200/220VAC 0.4kW 50/60 Hz 3A
<b>Outer Dimensions:</b>	250 x 600 x 770mm
<b>Instrument Weight:</b>	85kg



Sold & serviced by:

**OpTest Equipment Inc.**

[www.optest.com](http://www.optest.com) - [sales@optest.com](mailto:sales@optest.com) - +1-613-632-5169